

Accessories Cell Manipulation Calibre Scientific

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Cell Manipulation](#) / [Accessories Cell Manipulation Calibre Scientific](#)

Product Data

Item No.	description	Price
5192072001	Positioning aid, pack of 2, for mounting universal capillary holder on TransferMan® 4m/r and InjectMan® 4	CHF 128.00
5192080004	Y-cable FJ4, Y-cable for connecting FemtoJet 4i/4x with a PC and TransferMan 4r or InjectMan 4	CHF 158.00
5192081000	Y-cable PX, Y-cable for connecting PiezoXpert with TransferMan 4r or InjectMan 4	CHF 158.00
5196063007	Grip head set 4 (slim shape), for Capillary holder 4 slimshape, incl. 6 o-rings and 2 distance sleeves, for capillaries with outer diameter 1.0 mm	CHF 54.60
5196082001	Grip head set 4 size 0, incl. 6 o-rings and 2 distance sleeves, for capillaries with outer diameter 1.0 to 1.1 mm	CHF 64.68
5196083008	Grip head set 4 size 1, incl. 6 o-rings and 2 distance sleeves, for capillaries with outer diameter 1.2 to 1.3 mm	CHF 54.60
5196084004	Grip head set 4 size 2, incl. 6 o-rings and 2 distance sleeves, for capillaries with outer diameter 1.4 to 1.5 mm	CHF 54.60
5196085000	Grip head set 4 size 3, incl. 6 o-rings and 2 distance sleeves, for capillaries with outer diameter 0.7 to 0.9 mm	CHF 54.60
5196091000	Dust cap CellTram® 4, for filling valve (Oil-version)	CHF 45.00

5252070038	Connecting cable NI 2 Femto-Jet 4 with angled plug to connect a FemtoJet 4i/4x to an InjectMan NI 2	CHF 113.00
5252850130	Acoustic Silencer	CHF 108.00

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 20:04 o'clock.