

# Chemical Oxygen Demand (COD) Reagecon

[Produkt im eShop](#)

[Startseite](#) / [Marken](#) / [Reagecon Diagnostics](#) / [Chemical Oxygen Demand \(COD\) Reagecon](#)

## Produktdaten

Reagecon's offering includes a comprehensive range of COD Standards for use as Control Materials and COD Vials to be used in conjunction with Aqualytic® and Hach® spectrophotometers using the APHA Standard Method 5220D.

- **COD Standards** are ideal for use as Control Standards to verify that correct analysis for COD has taken place. Achieving an acceptable result for the Control Standard will improve confidence in sample readings for this important environmental parameter.
- **COD Reagents** include reagents for the two accepted methods for measuring COD at concentrations less than 400mg/l. Where the concentration is greater than 400mg/l, the sample must be diluted.
- **COD Reagent Vials** from Reagecon can be used in conjunction with the Aqualytic PC Spectro, PC Compact Vario and all Hach® spectrophotometers. This compatibility is proven in the Reagecon Technical Publication.

Art.-Nr.	Inhalt	Verpackung	Preis
500.420720	25 null	Glass Vials	CHF 100.70
500.420721	25 null	Glass Vials	CHF 100.70
500.420722	25 null	Glass Vials	CHF 100.70
500.AGS1H	500 Milliliter	Amber Glass Winchester Bottle	CHF 75.90
500.AGS1W	2 Liter	Amber Glass Winchester Bottle	CHF 158.30
500.COD10	500 Milliliter	HDPE Bottle GL45	CHF 195.80
500.COD100	500 Milliliter	HDPE Bottle GL45	CHF 155.40
500.COD1000	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.COD10M	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.COD1300	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.COD1500	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.COD20	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.COD200	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.COD2000	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.COD20M	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.COD3000	500 Milliliter	HDPE Bottle GL45	CHF 325.20
500.COD30K	1 Liter	HDPE Bottle GL45	CHF 210.55
500.COD50	500 Milliliter	HDPE Bottle GL45	CHF 131.55

500.COD500	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.COD5000	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.COD600	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.COD6000	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.COD60M5	500 Milliliter	HDPE Bottle GL45	CHF 131.55
500.CODMS	500 Milliliter	HDPE Bottle GL45	CHF 65.90
500.KC2002F	1 Liter	Amber Glass Winchester Bottle GL45	CHF 43.90
500.WTR50W	2 Liter	Amber Glass Winchester Bottle	CHF 121.90

\* Die Preise sind freibleibend und verstehen sich als Abgabepreise in Schweizer Franken ohne Mehrwertsteuer (MWST), sowie exklusive aller sonstigen Gebühren, Abgaben und Steuern. Die Preise, die im eShop angezeigt werden, sind aufgrund der regelmässigen Aktualisierung möglicherweise aktueller als auf der PDF-Datei.

\*\* Bitte beachten Sie, dass bei der Bestellung von Chemikalien und Reinigungsmittel Transport - und Verpackungskosten für Gefahrgüter sowie gesetzlich vorgeschriebene Gebühren berechnet werden. Diese werden Ihnen auf der Auftragsbestätigung, die Sie ergänzend zur Bestellbestätigung erhalten, detailliert ausgewiesen.

\*\*\* Weitere Informationen wie Technische Informationen und Sicherheitsdatenblätter finden Sie online in unserem eShop

\*\*\*\* Die PDF-Datei wurde erzeugt auf [www.huberlab.ch](http://www.huberlab.ch) am 18.05.2026 um 22:19 Uhr.