

## Accessories for Ultrasonic cleaners Elmasonic Elma



Product in eShop

[Home](#) / [Brands](#) / [elma](#) / Accessories for Ultrasonic cleaners Elmasonic Elma

### Product Data

Item No.	Description	UNSPSC	ECLASS	Price
626.0032.86	Elma Cover Typ 180 plastic, blue RAL 5002	41100000	32030000	CHF 45.00
626.0032.87	Elma Cover Typ 150/300 plastic, blue RAL 5002	41100000	32030000	CHF 55.00
626.0427.06	Insert basket modular size 300	41100000	32030000	CHF 1'780.00
626.0771.66	Elma insert Basket Typ 100 inox, 255 / 200 / 25 mm (Int.)	41100000	32030000	CHF 195.00
626.1004.24	Elma insert Basket with Rilsan-coated handles	41100000	32030000	CHF 190.00
626.1076.41	Stainless steel basket for Select900 with Rilsan-coated handles	41100000	32030000	CHF 465.00
626.1076.75	Basket for Elmasonic size 500 460 x 300 x 195mm, 11mm meshes	41100000	32030000	CHF 350.00
626.1089.16	Modular basket system made of stainless steel with handles made of plastic, for size 30	41100000	32030000	CHF 66.30
626.1118.38	Modular basket system made of stainless steel with handles made of plastic, for size 80	41100000	32030000	CHF 165.00
626.1119.01	Elma Mutebox ML, Sound Reductionbox Dim ext: 545 x 530 x 525 mm	41100000	32030000	CHF 1'090.00

626.1130.45	Modular basket system made of stainless steel with handles made of plastic, for size 40	41100000	32030000	CHF 75.00
626.1130.50	Elma Modular basket inox size 60, Dim: 255 / 105 / 95 mm	41100000	32030000	CHF 80.00
626.1130.72	Insert basket modular size 300 Dim. int: 450 / 255 / 135 mm	41100000	32030000	CHF 225.00
626.1133.87	Modular basket system made of stainless steel with handles made of plastic, for size 100	41100000	32030000	CHF 190.00
626.1133.90	Elma Modular basket inox size 120, Dim: 250 / 195 / 135 mm	41100000	32030000	CHF 175.00
626.1134.04	Modular basket system made of stainless steel with handles made of plastic, for size 180	41100000	32030000	CHF 232.20
626.1240.42	Elma Starter Kit for modular basket	41100000	32030000	CHF 60.00
626.1252.21	Elma Partition wall size S high for modular basket	41100000	32030000	CHF 46.00
626.2070.90	Elma acid-resistant insert tray for Elmasonic 100/120	41100000	32030000	CHF 180.00
626.2688.21	Elma Lid for Glas Ø90mm 600ml + 1000ml	41100000	32030000	CHF 9.85
626.8001.21	Elma tec clean A3, 10l-canister	41100000	32030000	CHF 115.00

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 15.05.2026 at 07:28 o'clock.