

Dräger X-act® 5000 Basic automatic tube pump



Product in eShop

Home / [Assortment](#) / [Safety & Hygiene](#) / Dräger X-act® 5000 Basic automatic tube pump

Product Data

The automatic tube pump Dräger X-act® 5000 Basic is the all-in-one solution designed for measurements with Dräger short-term tubes and for sampling tubes and systems. Ease of operation and a high degree of reliability compliment the measurement and sampling of gases, vapors and aerosols.

One order includes pump , battery and plug-in power supply.

Item No.	description	Price
18.4523.45	Plug-in power supply for Dräger X-act 5000 100-240VAC	CHF 98.90
18.4523.50	Dräger X-act 5000	CHF 1'962.00
18.4523.52	Accu NiMhy T 4 for Dräger X-act 5000	CHF 291.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 17:00 o'clock.