

## Snap Ring Vials ND11



Product in eShop

[Home](#) / [Assortment](#) / [Research & Development](#) / [Chromatography](#) / [Vials & flasks](#) / Snap Ring Vials ND11

### Product Data

#### The Universal Autosampler Vial

Universally usable vials for almost all autosamplers, even for those with robotic handling.

Vials can also be crimped with a standard 11mm Aluminum Crimp Seal, as the two snap ring lips have the same height as a crimp neck.

The vials are preferentially used on instruments of the following manufacturers: Agilent, CTC, DANI, Dionex, Jasco, Shimadzu, Spark, Thermo Scientific, Varian, VWR (Merck®)/Hitachi, Waters®, etc.

Item No.	description	vol.	h x d	colour	bottom	Price
----------	-------------	------	-------	--------	--------	-------

8.1102.00	Snap Ring Vial, 1st hydrol. class, wide opening	1,5 ml	32 x 11,6 mm	clear glass	flat	CHF 121.65
8.1102.00s	Snap Ring Vial, 1st hydrol. class, wide opening, silanized	1,5 ml	32 x 11,6 mm	clear glass	flat	CHF 460.90
8.1102.01s	Snap Ring Vial, 1st hydrol. class, wide opening, wide opening, label + filling lines, silanized	1,5 ml	32 x 11,6 mm	clear glass	flat	CHF 468.95
8.1102.02	Snap Ring Vial, 1st hydrol. class, wide opening, label + filling lines	1,5 ml	32 x 11,6 mm	amber glass	flat	CHF 138.70
8.1102.02s	Snap Ring Vial, 1st hydrol. class, wide opening, wide opening, label + filling lines, silanized	1,5 ml	32 x 11,6 mm	amber glass	flat	CHF 477.70
8.1102.03	Snap Ring Vial, 1st hydrol. class, wide opening, label + filling lines	1,5 ml	32 x 11,6 mm	clear glass	flat	CHF 130.20
8.1106.84	Snap Ring Vial with integrated 0,3 ml Mico-Insert, 1st hydrol class, Base Bonded	---	32 x 11,6 mm	clear glass	flat	CHF 734.85
8.1108.56	PP Snap Ring Micro-Vial	0,7 ml	32 x 11,6 mm	transparent	flat	CHF 198.05
8.1110.92	PP Snap Ring Micro-Vial	0,3 ml	32 x 11,6 mm	transparent	flat	CHF 198.05
8.1120.52	Snap Ring Vial with integrated 0.3 ml Micro-Insert, 1. hydrol. class, writing patch and filling lines, with widel micro-insert Base Bond	---	32 x 11.6 mm	glass clear	flat	CHF 771.40

8.1120.53	Snap Ring Vial with integrated 0.3 ml Micro-Insert, 1. hydrol. class, writing patch and filling lines, with widel micro-insert Base Bond	---	32 x 11.6 mm	amber glass	flat	CHF 795.10
-----------	--	-----	--------------	-------------	------	------------

## Accessories

Item No.	Description	Price
8.0602.18	Micro-Inserts for Short Thread Vials with 6 mm opening	from CHF 192.30
8.1101.68	PE Snap Ring Seals ND11	from CHF 69.25
8.1105.13	Septa for Snap Ring Seals ND11	from CHF 98.90

## More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 10.05.2026 at 10:33 oclock.

[www.huberlab.ch](http://www.huberlab.ch)