

Samplepack Headspace PerkinElmer

[Product in eShop](#)[Home](#) / [Brands](#) / [CG](#) / Samplepack Headspace PerkinElmer

Product Data

- 8.2002.38: 5 ml Headspace Vial , 38.2 x 22 mm, clear glass, 1st hydrol. class, HS-crimp neck, rounded bottom
- 8.2002.39: 20 ml Headspace Vial, 75.5 x 23 mm, clear glass, 1st hydrol. class, HS-crimp neck, rounded bottom
- 8.2002.49: 20 mm Headspace Cap with overpressure protection, clear lacquered, with hole, Pharma-Fix-Septa, Butyl/PTFE, 50° shore A, 3 mm
- 8.2006.41: 20 mm UltraClean Headspace Cap with overpressure protection, clear lacquered, with hole, silicone white/PTFE beige, 45° shore A, 3 mm
- 8.2002.44: 20 mm UltraClean Headspace Cap with overpressure protection, clear lacquered, with hole, silicone blue transparent/PTFE white, 45° shore A, 3 mm
- 8.2002.48: 20 mm Headspace Cap with overpressure protection, clear lacquered, with hole, moulded septa Butyl/PTFE grey, 50° shore A, 3 mm

Item No.	Price
8.2018.09S	for free

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 18.05.2026 at 22:22 o'clock.