

Snap Cap Vials ND18 + ND22 and Snap Caps



Product in eShop


[Home](#) / [Assortment](#) / [Research & Development](#) / [Chromatography](#) / [Vials & flasks](#) / Snap Cap Vials ND18 + ND22 and Snap Caps

Product Data

- Easy to handle and inexpensive storage vials with push-on Snap Caps
- Caps and Vials separately obtainable.
- Quickly and easily to reopen and reseal.
- No liners are required in the cap.
- For storage of powders and solids
- Appropriate PE Snap Caps ND18 + ND22, transparent, closed top

Item No.	description	h x d	colour	Price
8.1800.67	Snap Cap 18 mm PE , closed top, for ND18	5,2 x 19,8 mm	transparent	CHF 34.25
8.1801.89	Snap Cap Vial ND18	40 x 20 mm	clear glass	CHF 159.85
8.1801.90	Snap Cap Vial ND18	50 x 22 mm	clear glass	CHF 215.05
8.2200.68	Snap Cap 22 mm PE , closed top, for ND22	5,5 x 23 mm	transparent	CHF 35.40
8.2201.91	Snap Cap Vial ND22	48 x 26 mm	clear glass	CHF 262.90

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 09.05.2026 at 03:27 o'clock.