

Precision Balances PR Series Precision Ohaus

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Weighing](#) / [Scales](#) / Precision Balances PR Series Precision Ohaus

Product Data

Designed for Routine Weighing Applications in Your Workplace

Offering accuracy and repeatability in essential weighing applications in laboratory, industrial and education settings, PR Analytical and Precision Balances deliver competitive performance at an economical price. Featuring RS232 connectivity for easy communication, and a backlit display and a simple interface for uncomplicated operation, the PR is perfectly designed for your workplace.

Item No.	type no.	weighing range	display precision	pan	Price
16.4298.33	PR1602	up to 1600 g	0.01 g	180 mm	CHF 1'015.00
16.4298.34	---	---	---	---	CHF 1'170.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 21:12 o'clock.