

Power supply nanoPAC-300 Mini for gel electrophoresis tanks Clever Scientific

[Product in eShop](#)

Home / [Labortops Q1 2026](#) / Power supply nanoPAC-300 Mini for gel electrophoresis tanks Clever Scientific

Product Data

The nanoPAC-300 is capable of running all Clever Scientific horizontal multiSUB™ systems and vertical omniPAGE mini gel systems.

- Maximum constant current output of 400mA
- Simple two step feature to set a programmable change in voltage, current and time during the run
- Constant voltage up to 300 V
- Adjustable in 1 V and 1 mA increments
- Two parallel outputs
- Continuous run or timed setting up to 999 minutes
- Automatic safety shutdown
- Alarm function
- Easy-to-clean polycarbonate housing
- 3-digit LED display

Supply requirements 300 V / 400 mA

Item No.	type no.	Price
606.2727.11	nanoPAC-300P Mini	CHF 612.85

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 04.05.2026 at 04:23 o'clock.