

Volumetric trapezoidal flasks, borosilicate glass 3.3, class A LLG-Labware

[Product in eShop](#)

[Home](#) / [Labortops Q1 2026](#) / Volumetric trapezoidal flasks, borosilicate glass 3.3, class A LLG-Labware

Product Data

Trapezoidal shape. Conformity certified with batch ID.

- DIN EN ISO 1042 (600.6862.25, 600.6862.26, 600.6862.29, 600.6862.30, 600.6862.31)
- With standard ground joint and PE stopper
- Blue ring mark and inscription

Item No.	for nominal volume	Ground Joint	tol. +/-	Price
600.6862.25	1 ml	NS7/16	±0.025 ml	CHF 27.70
600.6862.26	2 ml	NS7/16	±0.025 ml	CHF 27.70
600.6862.27	3 ml	NS7/16	±0.025 ml	CHF 27.70
600.6862.28	4 ml	NS7/16	±0.025 ml	CHF 27.70
600.6862.29	10 ml	NS10/19	±0.040 ml	CHF 28.40
600.6862.30	20 ml	NS10/19	±0.040 ml	CHF 34.70
600.6862.31	25 ml	NS10/19	±0.040 ml	CHF 34.70

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 25.05.2026 at 15:35 o'clock.