



## Enduro™ power supplies

[Product in eShop](#)

Home / [Brands](#) / [Corning](#) / Enduro™ power supplies

### Product Data

Dependable performance, intelligent control, for the connection of different chamber systems. Ideal for DNA, RNA and protein electrophoresis as well as blotting.

- Simple to program
- Select constant voltage or constant current mode
- Compact, lightweight and stackable
- Extremely quiet operation

#### Specifications

Output voltage: 2 ... 300 V/1 V

range/increments: 4 ... 500 mA/1 mA 90 W max/na

Timer: 1 min to 99 h 59 min or continuous

Programmable: No

Dimensions: 190 x 250 x 80 mm

Weight: 2.2 kg

Input voltage: 120-230 V

Item No.	description	Price
606.5953.14	Enduro 300 V power supply	CHF 692.00

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 21.05.2026 at 00:55 o'clock.