

## Capillary columns GC BPX90 Trajan

[Product in eShop](#)

Home / [Labortops Q2 2026](#) / Capillary columns GC BPX90 Trajan

### Product Data

90% Cyanopropyl Polysilphenylene-siloxane. For the fast separation of fragrances, aromatics, pesticides, PCB, isomers of fatty acid methyl esters.

- Bonded phase
- Highly polar
- Thermally stable
- Able to be solvent rinsed

Temp. range: 80 ... 260 / 280 °C

Item No.	Description	UNSPSC	ECLASS	Price
619.2251.62	Capillary column BPX90 60 m x 0,25 mm ID gauge 0,25 µm	41100000	32030000	CHF 1'050.00 <del>CHF 2'290.00</del>
619.6824.18	Capillary column GC BPX90 0,25mm ID x 0,25µm x 15m length	41100000	32030000	CHF 642.00
619.6824.19	Capillary column GC BPX90 0,25mm ID x 0,25µm x 30m length	41100000	32030000	CHF 1'015.00
619.6824.20	Capillary column GC BPX90 0,32mm ID x 0,5µm x 30m length	41100000	32030000	CHF 1'050.00

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

[www.huberlab.ch](http://www.huberlab.ch)