

Biohazard tabletop waste container neoLab®



Product in eShop

[Home](#) / [Assortment](#) / [Safety & Hygiene](#) / [Disposal](#) / [Disposal containers](#) / Biohazard tabletop waste container neoLab®

Product Data

Automatic hands-free operation for the disposal of biohazardous waste

The compact table-top design is ideal for odour control and is perfect for small items such as pipette tips and microcentrifuge tubes.

The automatic lid opening mechanism makes the unit convenient for hands-free use and closes automatically and securely to ensure a more hygienic environment.

Features include a polypropylene inner liner with a handle for easy emptying of the contents.

The stainless steel body and chrome-plated lid are not only aesthetically pleasing but also durable. They are resistant to bleach and easy to clean.

The appliance requires three AAA batteries, which are not included.

CE and RoHS certification are available.

- 1.5 litre volume
- 17 cm high
- 12.4 cm top diameter

Item No.	description	Price
19.6215.19	Biohazard tabletop waste container	CHF 77.35

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 06:37 o'clock.

www.huberlab.ch