

# Saveris 2-T2 testo wireless temperature logger



 Product in eShop

Home / [Brands](#) / [testo](#) / Saveris 2-T2 testo wireless temperature logger

## Product Data

The WiFi data logger is ideal for monitoring temperature in rooms. Measurement data is sent to the testo Cloud (online data storage) via wireless LAN and can be accessed and evaluated on all internet-enabled end devices. The testo Cloud Basic package can be accessed free of charge.

- Internal data memory up to 10000 measured values per channel
- Limit values adjustable via testo-Cloud
- E-mail notification when limit values are exceeded
- Current measured values can be read on the display
- Micro USB interface
- 2 connections for external NTC temperature probes

**At least one NTC temperature probe is required to operate the testo Saveris 2-T2. Please order separately.**

**Scope of supply:** WiFi Temperature logger, testo-Cloud access (basic functionality), micro USB cable, wall bracket, batteries, calibration protocol

- Sensor type: NTC
- Temperature range: -50 ... 150 °C
- Accuracy:  $\pm 0.3$  °C
- Measuring cycle: 15 min ... 24 h
- Communication cycle: 15 min ... 24 h
- Memory: 10000 measured values per channel
- Dimensions: 95 x 75 x 30.5 mm
- Weight: 240 g
- Power supply: 4 x AA AIMn Mignon batteries
- IP code: IP65
- Measuring range temperature: 150 °C
- Storage readings: 10000

Item No.	type no.	Price
19.5722.32	testo Saveris 2-T2	CHF 250.00

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 18.05.2026 at 22:21 o'clock.

[www.huberlab.ch](http://www.huberlab.ch)