

Luminometer Clean-Trace™ LM1 3M

[Product in eShop](#)[Home](#) / [Brands](#) / [3M](#) / Luminometer Clean-Trace™ LM1 3M

Product Data

Simple, quick and reliable hand held diagnostic tool to verify cleaning effectiveness and capture, store and manage data for later analysis.

- Compact and easy-to-use
- Rapid results allow immediate corrective action
- To be used with the dry swabs 3M™ Clean-Trace™ UXL100, AQT200 and AQF100
- Do not use outdoors

- Operating temperature range: +5 °C to +40 °C (+41 °F to +104 °F)
- Relative humidity: 10 - 90 %, non-condensing
- Storage temperature range: -40 °C to +70 °C (-40 °F to +158 °F)
- AC voltage for power supply unit: 100 - 240 V
- External connections: USB 2.0 type mini AB connector
- Width: 84.47 mm
- Height: 268.17 mm
- Weight: 517.83 g

Item No.	Price
17.9995.10	CHF 4'202.00

Accessories

Item No.	Description	Price
12.9995.02	Dry swabs for Luminometer 3M™ Clean-Trace™ NG3/LM1 3M	from CHF 495.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.05.2026 at 03:50 oclock.

www.huberlab.ch