

Greenough- Stereomicroscope Educational-Line OSE Kern

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Quality Assurance](#) / [Microbiology](#) / [Greenough-Stereomicroscope Educational-Line OSE Kern](#)

Product Data

Ideal for workshops, schools and training.

- Tube 45° inclined
- WF 10x/ Ø 20 mm
- Diopter adjustment: Both-sided
- Illumination 1 W LED incident and transmitted
- Integrated handle and arm curved stand

Total dimension (W x D x H): 230 x 130 x 300 mm

Weight: 2.0 kg

Power supply: 230 V / 50 Hz, battery

Magnification: 2x / 4x

Item No.	type no.	description	Price
17.4515.40	OSE 421	Binocular	CHF 391.50

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 00:19 o'clock.