

Exhibition glass Duran® DWK

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Quality Assurance](#) / [Sample preparation](#) / [Glassware](#) / [Exhibition glass Duran® DWK](#)

Product Data

- Manufactured from borosilicate glass 3.3
- Overhead position of the specimen jar
- Tighter closure due to precise grinding
- Fired-on, highly durable white ceramic

The multitude of different work steps in the laboratory environment demands a multitude of universally applicable glassware. These glass devices usually have a very simple shape and are characterized by their versatility. Exhibition glasses and museum jars are used to store materials and are often used in exhibitions. Here the objects on display can be preserved and then viewed without touching them.

Item No.	capacity	diameter	height	Price
8.8024.01	100 ml	52 mm	135 mm	CHF 917.00
8.8039.01	300 ml	69 mm	163 mm	CHF 1'223.00
8.8048.01	600 ml	81 mm	214 mm	CHF 1'454.50
8.8051.01	750 ml	90 mm	240 mm	CHF 1'587.50
8.8058.01	1200 ml	100 mm	253 mm	CHF 2'087.50

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 04.05.2026 at 13:36 o'clock.