

Cylinder with knobbed lid Duran® DWK

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Quality Assurance](#) / [Sample preparation](#) / [Glassware](#) / [Cylinder with knobbed lid Duran® DWK](#)

Product Data

- Manufactured from borosilicate glass 3.3
- With polished rim
- The lid can be opened easily while using the button
- Fired-on, highly durable white ceramic

Ideally, containers that can be closed are used to store substances and materials. In this way, substances and materials can be protected from unwanted effects.

| Item No. | capacity | diameter | height | Price |
|-----------|----------|----------|--------|--------------|
| 8.5010.01 | 250 ml | 80 mm | 80 mm | CHF 853.00 |
| 8.5030.01 | 500 ml | 100 mm | 100 mm | CHF 1'074.50 |
| 8.5050.01 | 1000 ml | 120 mm | 120 mm | CHF 114.85 |
| 8.5100.01 | 2000 ml | 150 mm | 150 mm | CHF 157.50 |
| 8.5210.01 | 6000 ml | 210 mm | 210 mm | CHF 305.15 |
| 8.5320.01 | 12000 ml | 260 mm | 260 mm | CHF 507.75 |

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 02.05.2026 at 04:41 o'clock.