

# Glass Tubes for NMR analysis TEOtube Zeochem



 Product in eShop

Home / [Brands](#) / [zeochem](#) / Glass Tubes for NMR analysis TEOtube Zeochem

## Product Data

### Glass Tubes for NMR analysis

- Made of Borosilicate Glass 5.1 ASTM Type 1 class B glass
- Tubes with orange, sky blue or black caps
- Length: 177.8 mm
- Wall thickness 0.43 mm

### Tube with orange cap

- TA - Throw Away
- Outer Diameter  $4.950 \pm 0.030$  mm
- Concentricity: 25 m
- Camber: 75 m

### Tube with sky blue cap

- HP - High Precision
- Outer Diameter:  $4.950 \pm 0.015$  mm
- Concentricity 7 m
- Camber 13 m

### Tube with black cap

- UP - Ultra Precision
- Outer Diameter:  $4.950 \pm 0.010$  mm
- Concentricity 4 m
- Camber 4 m

Item No.	type no.	Price
198.4000.42	Cap orange	CHF 16.30
198.4000.51	Cap sky blue	CHF 16.30

## More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 13.05.2026 at 06:32 o'clock.