

PFAS Safe Guard 2in1 Kits - Snap ring vials and closures



Product in eShop

[Home](#) / [Assortment](#) / [Research & Development](#) / [Chromatography](#) / [Vials & flasks](#) / [PFAS Safe Guard 2in1 Kits - Snap ring vials and closures](#)

Product Data

PFAS stands for per- and polyfluoroalkyl substances and is colloquially known as 'forever chemicals'. They are a group of man-made chemicals that have been used in various industrial and consumer products since the 1940s. PFAS are known for their strong chemical bonds that make them resistant to heat, water and oil. This property has led to their use in a variety of applications, including firefighting foams, non-stick cookware (such as Teflon), waterproof fabrics, food packaging and many other household and industrial products.

PFAS can be found in:

- Environmental monitoring
- product testing
- Industrial quality assurance
- Commercial household cleaning products
- Food packaging and much more.

The EPA's goal is to ensure science-based decision making to identify which PFAS may pose risks to human health and the environment and to develop tests to measure, remove and destroy them.

Using Chroma Globe PFAS sample analyser caps ensures PFAS-free results!

When taking PFAS samples, it is important that vials and caps are PFAS-free. Using Chroma Globe's new PFAS caps ensures the integrity of the data collected and fulfils all required collection standards. Chroma Globe has developed a range of caps that are compatible with all standard plastic screw-top or snap-top bottles. These caps are recommended for the testing of drinking water, food packaging and commercial household cleaning products.

Item No.	description	packaging	Price
----------	-------------	-----------	-------

8.1122.67	0.3ml PP snap ring micro vial, 32 x 11.6mm, transparent (8.1110.92) + PFAS Safe Guard PE snap ring cap, blue, silicone white/polypropylene transparent	100 vials and 100 caps in a white 2in1 kit box	CHF 50.60
8.1122.68	0.3ml PP snap ring micro vial, 32 x 11.6mm, transparent (8.1110.92) + PFAS Safe Guard PE snap ring cap, blue, silicone white/polypropylene transparent, slotted	100 vials and 100 caps in a white 2in1 kit box	CHF 52.90
8.1122.69	0.7ml PP snap ring micro vial, 32 x 11.6mm, transparent (8.1108.56) + PFAS Safe Guard PE snap ring cap, blue, silicone white/polypropylene transparent	100 vials and 100 caps in a white 2in1 kit box	CHF 50.60
8.1122.70	0.7ml PP snap ring micro vial, 32 x 11.6mm, transparent (8.1108.56) + PFAS Safe Guard PE snap ring cap, blue, silicone white/polypropylene transparent, slotted	100 vials and 100 caps in a white 2in1 kit box	CHF 52.90

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 29.04.2026 at 22:43 o'clock.

www.huberlab.ch