

Sample bags, TWIRL'EM, without writing area, extra heavy sample closure Labplas



Product in eShop

Home / Assortment / Medicine & Healthcare / Sampling / Sampling bags / Sample bags, TWIRL'EM, without writing area, extra heavy sample closure Labplas

Product Data

For collecting, transporting and analysing liquid and solid samples.

- Integrated secure closure system with 2 flat wires
- Sterile
- Contaminant-free
- Waterproof and airtight
- RNase, DNase and pyrogen free
- Without labelling field

Item No.	volume	Dimensions	Thickness	Price
15.1123.80	5500 ml	305 x 406 mm	0.101 mm	CHF 274.00
15.1123.81	6750 ml	305 x 457 mm	0.101 mm	CHF 300.00
15.1123.82	8000 ml	305 x 508 mm	0.101 mm	CHF 395.00
15.1123.83	10500 ml	305 x 610 mm	0.101 mm	CHF 448.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 25.04.2026 at 12:24 o'clock.