

Mini Diaphragm Vacuum Pump LABOPORT® N 938.50 KT.18 KNF

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Liquid Handling](#) / [Pumps](#) / [Mini Diaphragm Vacuum Pump LABOPORT® N 938.50 KT.18 KNF](#)

Product Data

Chemically resistant, oil-free diaphragm vacuum pump used for a wide range of laboratory applications. It features PPS heads, PTFE-coated diaphragms and FFPM valves

- Flow Rate (max.): 30 l/min
- Ultimate Vacuum (max.): 15 mbar (abs.)
- Clean, 100% oil-free operation
- Corrosion-resistant
- Compact, quiet and reliable
- Small footprint

Item No.	type no.	Price
17.8621.38	N 938.50 KT.18	CHF 2'086.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 20:10 o'clock.