


Printer dataLink® pro interscience

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Life Science](#) / [Cell Biology](#) / [Colony counters](#) / [Printer dataLink® pro interscience](#)

Product Data

dataLink® pro is a traceability system for the microbiological analysis workflow. It allows the transmission, printing and labeling on the edge of the Petri dish of the sample. dataLink® pro is connected to the easySpiral® platers and the Scan® automatic colony counters.

- Semi-automated labelling
- Traceability of data
- Easy counting
- Fast labelling
- Three label formats: sample number, barcode, datamatrix
- Label size: 5 x 50 mm

Scope of supply:

Printer, power cord, USB cable, 2 user manuals, thermal transfer ribbon, mandrel, 3 x 2500 adhesive labels, USB with the monitoring software

- Weight: 12 kg
- Dimensions (WxDxH): 252 x 288 x 460 mm
- Rated capacity: 100 W
- Power supply: 100 ... 240 V, 50 ... 60 Hz
- Warranty: 1 year

| Item No. | type no. | Price |
|------------|---------------|---------------|
| 17.4370.10 | dataLink® pro | CHF 11'882.20 |

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 02.05.2026 at 08:17 oclock.

www.huberlab.ch