

Short thread vials ND9, PP LLG-Labware

[Product in eShop](#)

[Home](#) / [Labortops Q4 2025](#) / Short thread vials ND9, PP LLG-Labware

Product Data

The LLG-Short thread vials are suitable for PFAS analyses.

- Made of PP
- Free from polyhaloolefins

Colour: Transparent

With labelling field: No

Item No.	vol.	h x d	execution	Price
600.0082.62	0.3 ml	32 mm x 11.6 mm	Inner cone, 0.3 ml	CHF 21.30
600.7125.14	1.5 ml	32 mm x 11.6 mm	with filling mark	CHF 22.60
600.7125.15	0.7 ml	32 mm x 11.6 mm	Micro-insert 0.7 ml	CHF 21.30

Accessories

Item No.	Description	Price
600.7125.16	Short thread seals ND9, PP, ready assembled LLG-Labware	from CHF 30.60

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 16:59 o'clock.

www.huberlab.ch