

## Laboratory beaker, stainless steel Bürkle

[Product in eShop](#)

[Home](#) / [Labortops Q1 2026](#) / [Laboratory beaker, stainless steel Bürkle](#)

### Product Data

The laboratory beaker is suitable for sterile environments and food technology applications. Its universal Griffin beaker shape makes it the ideal choice for stirring, mixing or sample preparation.

- With spout and rim
- Smooth surface, easy to clean, hygienic
- Contact base for optimum heat distribution
- Corrosion and chemical resistant
- Unbreakable
- Autoclavable
- Food-safe

Item No.	volume	diameter	height	Price
15.7121.01	100 ml	50 mm	60 mm	CHF 24.80
15.7121.02	250 ml	65 mm	85 mm	CHF 27.10
15.7121.03	500 ml	80 mm	110 mm	CHF 42.90
15.7121.04	1000 ml	100 mm	130 mm	CHF 63.30
15.7121.05	3000 ml	150 mm	180 mm	CHF 94.90
15.7121.06	5000 ml	180 mm	210 mm	CHF 101.65

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 04.05.2026 at 08:28 oclock.

[www.huberlab.ch](http://www.huberlab.ch)