

## Sound proof box LS 40 Bandelin

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Microbiology](#) / [Homogenizers](#) / [Sound proof box LS 40 Bandelin](#)

### Product Data

To reduce the noise level when sonicating samples with an ultrasonic homogeniser. The lower part of the interior and the removable bottom tray are made of stainless steel and are therefore chemically resistant and easy to clean.

- Noise reduction of up to 30 dB-AU
- LED interior lighting
- Acrylic glass window
- With effective ventilation system
- With opening (40 x 70 mm) for cooling tubing, sample circulation or the SONOPULS temperature sensor on the back
- Removable bottom tray
- Easy retrofitting for use with sonotrodes or cup booster possible

**The stand HG 40 (article number: 17.4400.05) is required to attach the ultrasonic converter. Please order separately.**

- Int. Dim. (WxDxH): 350 x 310 x 880 mm
- Power supply: 100 ... 240 V, 50/60 Hz

Item No.	type no.	Price
17.3682.20	LS 40 for HD 2000.2 / HD 3000 / HD 4000	CHF 2'595.00

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 18.05.2026 at 21:09 oclock.

[www.huberlab.ch](http://www.huberlab.ch)