

Dispenser box KEK

[Product in eShop](#)

[Home](#) / [Labortops Q2 2026](#) / [Dispenser box KEK](#)

Product Data

- Material stainless steel (1.4301 / AISI 304)
- Functional surfaces in KEK-Finish
- To dispense ofdisposable items (hoods, gloves, etc.)
- 1 - 4 boxes
- Hinged lid with bow handle for filling
- Viewing slots on the front to monitor the filling level
- Embossed removal opening
- Removable drip tray arranged under the removal opening
- 2 mounting holes on the back for wall mounting

Item No.	type no.	description	dimensions (w x d x h)	Price
6105.3284.20	SPB-2-CLN	No. of compartments: 2	300 x 200 x 600 mm	CHF 355.00 CHF 630.00
6105.4949.13	SPB-4-CLN	No. of compartments: 4	600 x 200 x 600 mm	CHF 490.00 CHF 890.00
6105.4949.14	SPB-1-CLN	No. of compartments: 1	150 x 200 x 600 mm	CHF 285.00 CHF 490.00
6105.4949.15	SPB-3-CLN	No. of compartments: 3	450 x 200 x 600 mm	CHF 425.00 CHF 750.00

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 24.04.2026 at 01:47 o'clock.

www.huberlab.ch