

# Mobile Phase Reagecon

[Product in eShop](#)[Home](#) / [Brands](#) / [Reagecon Diagnostics](#) / [Mobile Phase Reagecon](#)

## Product Data

Reagecon's range of pre-mixed ready to use Mobile Phase Solutions are designed to be used as Buffered Eluents for Liquid Chromatography. The use of these high quality buffer systems will minimise variations of mobile phase pH, leading to dramatically improved selectivity, retention factor, peak shape, resolution and reproducibility.

Item No.	Content	Packaging	Price
500.MPFORM01W	2 liter	Amber Glass Winchester Bottle	CHF 100.80
500.MPFORMETH	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 42.15
500.PCMP7-W	2 liter	Amber Glass Winchester Bottle	CHF 234.80
500.SFABUFW	2 liter		CHF 44.85
500.SFAMIXW	2 liter		CHF 54.47

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 04.05.2026 at 08:28 o'clock.