

Petrochemical Gas Calibration Standards Reagecon

[Product in eShop](#)

Home / [Reagecon PetrochemicalStandards](#) / Petrochemical Gas Calibration Standards Reagecon

Product Data

Reagecon offers a range of Petrochemical Standards. These ready to use standards are used in fuel analysis as either a calibrant or analytical control solution. They can also be used to validate appropriate test methods or qualify an instrument in a regulated industry.

Item No.	Content	Packaging	Price
500.REGASCAL-1-250	250 milliliter	Amber Glass Winchester Bottle	CHF 211.35
500.REGASCAL-2-250	250 milliliter	Amber Glass Winchester Bottle	CHF 1'592.40
500.REGASCAL-5-250	250 milliliter	Amber Glass Winchester Bottle	CHF 704.15
500.REGASCAL-6-250	250 milliliter	Amber Glass Winchester Bottle	CHF 702.95
500.REGASCAL-7-250	250 milliliter	Amber Glass Winchester Bottle	CHF 702.95
500.REGASCAL-8-250	250 milliliter	Amber Glass Winchester Bottle	CHF 702.95

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

www.huberlab.ch