

Conductivity Standards ISO 17034 Certified Reference Materials (CRMs) Reagecon



Product in eShop

[Home](#) / [Brands](#) / [Reagecon Diagnostics](#) / [Conductivity Standards ISO 17034 Certified Reference Materials \(CRMs\) Reagecon](#)

Product Data

Reagecon manufacture a range of Conductivity Certified Reference Materials (CRMs), which are manufactured and certified in accordance with the requirements of ISO 17034. Reagecon's Conductivity Standards are all aqueous based, thereby eliminating any errors attributable to matrix mismatch. Standards have a specification of $\pm 1\%$ this high level of accuracy enables the Standards to be used as calibrators and/or controls in fulfilment of the most exacting industrial statutory requirements.

Item No.	Description	UNSPSC	ECLASS	Price
500.CRMCSKC100	Conductivity Standard 100 $\mu\text{S}/\text{cm}$ at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 255.00
500.CRMCSKC1000	Conductivity Standard 1000 $\mu\text{S}/\text{cm}$ at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 255.00
500.CRMCSKC10M	Conductivity Standard 10,000 $\mu\text{S}/\text{cm}$ at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 255.00
500.CRMCSKC20	Conductivity Standard 20 $\mu\text{S}/\text{cm}$ at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 255.00

500.CRMCSKC20 0	Conductivity Standard 200 µS/cm at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 255.00
500.CRMCSKC2 M	Conductivity Standard 2000 µS/cm at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 255.00
500.CRMCSKC3 M	Conductivity Standard 3000 µS/cm at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 255.00
500.CRMCSKC50	Conductivity Standard 50 µS/cm at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 255.00
500.CRMCSKC50 0	Conductivity Standard 500 µS/cm at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 255.00
500.CRMCSKC5 M	Conductivity Standard 5000 µS/cm at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 255.00
500.CRMCSKC84	Conductivity Standard 84 µS/cm at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 64.80
500.CRMCSKCL	Conductivity Standard 1413 µS/cm at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 64.80
500.CRMCSKCS	Conductivity Standard 147 µS/cm at 25°C ISO 17034 Certified Reference Material 500ml	41100000	32030000	CHF 64.80

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 11.05.2026 at 12:07 oclock.

www.huberlab.ch