

Technical pH buffer solutions, certified reference material (CRM) Reagecon



 Product in eShop

[Home](#) / [Brands](#) / [Reagecon Diagnostics](#) / Technical pH buffer solutions, certified reference material (CRM) Reagecon

Product Data

Certified reference material according to ISO 17034. The pH buffer solutions are colour-coded and the uncertainty of measurement of ± 0.02 pH was tested at 25 °C. Supplied with detailed certification report which includes traceability to NIST. Lot numbers and expiry dates are printed on the label.

- Capacity: 500 ml

Item No.	Content	Packaging	Price
500.CRMTB100	1 piece		CHF 37.05
500.CRMTB401	1 piece		CHF 37.05
500.CRMTB460	1 piece		CHF 37.05
500.CRMTB700	1 piece		CHF 37.05
500.CRMTB921	1 piece		CHF 37.05

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 03:17 o'clock.

www.huberlab.ch