

Osmolality Standards Reagecon

[Product in eShop](#)[Home](#) / [Labortops Q2 2026](#) / [Osmolality Standards Reagecon](#)

Product Data

Reagecon's range of Osmolality Standards for Osmometry determinations start at 50 mOsm/Kg and go up to 3000 mOsm/Kg. These standards are used for calibration, quality control, method validation and instrument calibration. They are manufactured in accordance with European and United States Pharmacopoeia guidelines. Products specifically designed for clinical applications are also available.

Item No.	Content	Packaging	Price
500.RE-OSM-0	1 piece	12x5ml Amber Glass Ampoule	CHF 107.00 CHF 141.10
500.RE-OSM-100	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM-1000	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM-1100	1 piece	12x5ml Amber Glass Ampoule	CHF 107.00 CHF 141.10
500.RE-OSM-1500	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM-200	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM-2000	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM-2500	1 piece	12x5ml Amber Glass Ampoule	CHF 207.00 CHF 274.15
500.RE-OSM-290	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM-300	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM-3000	1 piece	12x5ml Amber Glass Ampoule	CHF 207.00 CHF 274.15
500.RE-OSM-3500	1 piece	12x5ml Amber Glass Ampoule	CHF 110.00 CHF 145.30
500.RE-OSM-400	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM-50	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM-500	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM-600	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20

500.RE-OSM-700	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM-850	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM-900	1 piece	12x5ml Amber Glass Ampoule	CHF 113.00 CHF 148.20
500.RE-OSM10001	100 milliliter	HDPE Bottle	CHF 220.00 CHF 289.55
500.RE-OSM50001	100 milliliter	HDPE Bottle	CHF 220.00 CHF 289.55
500.RE-POSM-240	1 piece	12x5ml Amber Glass Ampoule	CHF 148.20
500.RE-POSM-280	1 piece	12x5ml Amber Glass Ampoule	CHF 274.15
500.RE-POSM-320	1 piece	12x5ml Amber Glass Ampoule	CHF 148.20
500.RE-ROSM-300	1 piece	12x5ml Amber Glass Ampoule	CHF 148.20
500.RE-ROSM-800	1 piece	12x5ml Amber Glass Ampoule	CHF 148.20

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 02.05.2026 at 04:43 o'clock.