

# Accessories for Combitips advanced® Eppendorf

[Product in eShop](#)

Home / [Labortops Q1 2026](#) / Accessories for Combitips advanced® Eppendorf

## Product Data

Adapter colour suits the colour of Combitips®. Marking of the adapter shows the size of Combitips®.

### Optimised **adapter connection:**

- Good optical feedback: 'window' shows correct closure
- Good haptic feedback: tight fit between adapter and Combitip®

The **Combitip®-rack** features trouble-free and single-hand operation during dispensing work. The transparent tray allows for the sizes of the Combitip® to be easily identified using the colour coding.

- Resealable, reusable box, autoclavable (Combitips® are not autoclavable)

Item No.	execution	designation	colour	Price
11.4850.88	Biopur®	adapter, 25 ml (individually wrapped)	red	CHF 103.25
11.4850.89	Eppendorf Quality	Adapter, 25 ml	red	CHF 14.15
11.4850.90	Eppendorf Quality	Adapter, 50 ml	light grey	CHF 14.15
11.4850.91	Biopur®	Adapter, 50 ml	light grey	CHF 103.25

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 03.05.2026 at 06:10 o'clock.