

Taschen-Thermometer, rote Spezial- füllung, Messb. -35° bis +50°C :1.0

[Product in eShop](#)[Home](#) / [Deals and news](#) / [Sale](#) /

Product Data

- ca. 140x11,5mm in vernickelter Doppelsteckhülse
- Einschlussform, ganz eintauchend
- Füllung: Quecksilber / rote Flüssigkeit / blaue Flüssigkeit

Komplett : Thermometer und Hülse
Ersatzthermometer **einzel**n

Item No.	execution	measuring range °C	length	filling	division	Price
4.4034.02	komplett	-35 +50 °C	145 mm	rot	1,0 °C	CHF 33.44
4.4034.04	komplett	-10 +100 °C	145 mm	rot	1,0 °C	CHF 15.00 CHF 34.31
4.4034.06	komplett	-10 +250 °C	145 mm	blau	2,0 °C	CHF 39.86
4.4036.02	einzel	-35 +50 °C	140 mm	rot	1,0 °C	CHF 18.16
4.4036.04	einzel	-10 +100 °C	140 mm	rot	1,0 °C	CHF 17.30
4.4036.06	einzel	-10 +250 °C	140 mm	blau	2,0 °C	CHF 29.28
4.4038.02	komplett	-35 +50 °C	145 mm	Hg	1,0 °C	CHF 35.65
4.4038.04	komplett	-10 +100 °C	145 mm	Hg	1,0 °C	CHF 29.56
4.4038.08	komplett	-10 +250 °C	145 mm	Hg	2,0 °C	CHF 36.23
4.4040.02	einzel	-35 +50 °C	140 mm	Hg	1,0 °C	CHF 20.02
4.4040.04	einzel	-10 +100 °C	140 mm	Hg	1,0 °C	CHF 19.18
4.4040.08	einzel	-10 +250 °C	140 mm	Hg	2,0 °C	CHF 25.68

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

www.huberlab.ch