

Schutzbeutel 250mm lang, 150mm breit, Stärke 50µ mit flüssigkeitsdichtem Sicherheitsverschluss



Product in eShop

[Home](#) / [Assortment](#) / [Research & Development](#) / [Microbiology](#) / [Bags](#) /

Product Data

Sekundärverpackung für die Beförderung von diagnostischen Proben entsprechend der Gefahrgutverordnung.

- Aus extrem reissfester Polyethylen-Spezialfolie mit vorfixiertem Sicherheitsverschluss, für diagnostische Proben der Risikogruppen 1 bis 3
- Sie sind flüssigkeitsdicht (BAM-Zulassung) im Sinne der Gesetzgebung (ADR, Verpackungsverordnung P 620/650) und verhindern Infektionsgefahren
- Mit fortlaufender Nummerierung, Barcode und eindeutiger Kennzeichnung
- Mit 30 mm Sicherheitsverschluss, Stärke 70 µm

Item No.	external dimensions	inner dimensions	execution	Price
3.6449.01	250 x 150 mm	240 x 120 mm	-	CHF 43.17
3.6449.02	300 x 200 mm	290 x 170 mm	-	CHF 67.86
3.6449.03	250 x 150 mm	240 x 120 mm	inkl. Saugelinlage	CHF 50.05
3.6449.04	300 x 200 mm	290 x 170 mm	inkl. Saugelinlage	CHF 101.00

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 08.05.2026 at 11:44 oclock.

www.huberlab.ch