

FILTER-BIO PES Sterile syringe filters

[Product in eShop](#)

Home / [FILTER-BIO](#) / FILTER-BIO PES Sterile syringe filters

Product Data

- in Polyestersulfone (PES)
- sterilized by gamma radiation
- individually packed
- every Filter is printed with expiry date

Item No.	pore size	diameter	Price
85.2000.00	0.22 µm	13 mm	CHF 135.65
85.2000.01	0.45 µm	13 mm	CHF 135.65
85.2000.02	0.22 µm	25 mm	CHF 80.90
85.2000.03	0.45 µm	25 mm	CHF 80.90
85.2000.04	0.22 µm	30 mm	CHF 95.20
85.2000.05	0.45 µm	30 mm	CHF 95.20

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 21:19 o'clock.