

Labcoat with press studs Tyvek® 500 model PL30NP DuPont™

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Safety & Hygiene](#) / [Personal Safety](#) / [Body protection](#) / Labcoat with press studs Tyvek® 500 model PL30NP DuPont™

Product Data

- Certified according to Regulation (EU) 2016/425
- Partial body chemical protective clothing, Category III, Type PB [6]
- EN 14126 (barrier to infective agents)
- Antistatic treatment (EN 1149-1) - on both sides

Item No.	size	Price
18.0411.01	M	CHF 361.50
18.0411.02	L	CHF 361.50
18.0411.03	XL	CHF 361.50
18.0411.04	XXL	CHF 361.50

Accessories

Item No.	Description	Price
18.0414.02	Sleeve Tyvek® 500 model PS32LA DuPont™	from CHF 208.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 18.05.2026 at 22:21 o'clock.

www.huberlab.ch