

neodisher™ LM 3



Product in eShop

[Home](#) / [Assortment](#) / [Chemicals](#) / [Cleaning Chemicals](#) / [Manual cleaning](#) / neodisher™ LM 3

Product Data

Alkaline cleaner for laboratory glass, liquid concentrate.

**For manual cleaning and pre-cleaning of laboratory glassware.
For the cleaning of pipettes and viscometers in special washers.**

Especially suitable for the cleaning of glassware in water, phosphates and metal laboratories and laboratories where enzymatic and serological analyses are conducted.

Also suitable for automated cleaning.

- Free of phosphates
- Low foam
- UN 1814 KALIUMHYDROXIDLÖSUNG, 8, II, E

Item No.	execution	capacity	Price
18.9410.01	bottle	12 x 1 liter	CHF 324.60
18.9410.02	canister	10 liter	CHF 216.35
18.9410.03	bottle	1 liter	CHF 27.20

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.05.2026 at 05:16 oclock.

www.huberlab.ch