

EasyDip™ Slide Staining System

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Medicine & Healthcare](#) / [Microscopy](#) / [Slides](#) / EasyDip™ Slide Staining System

Product Data

- Made of acetal
- Not autoclavable. Not for use with phenol, iodine, or ferric chloride. Not recommended for Gram stain.

EasyDip™ slide staining rack

The EasyDip™ Slide Staining Rack will hold up to 12 microscope slides with dimensions such as 75 x 25 mm (3 x 1 in.) or 76 x 26 mm with a thickness of 1.0 and 1.2 mm. The slides fit into individual slots for free passage and rapid drainage of staining fluids. Available in dark gray only. EasyDip™ jars sold separately.

- Holds up to 12 microscope slides
- Individual slide slots
- Easy slide removal

EasyDip™ Slide staining jars

Made of acetal, the EasyDip™ Staining Jar will not break like most glass jars do. It will resist most staining agents including alcohol and xylene (but not phenol, iodine, or ferric chloride). The wide stable base offers greater stability while the inside is recessed, allowing for a smaller reagent volume of only 80 ml.

Easy to clean and no metals to corrode. Ideal for special stains, frozen sections, and special processes.

Available in 5 different colors to help better identify contents or applications. Racks are sold separately.

- Will resist most commonly used histological solvents
- 80 ml reagent volume
- Individual jars are available in 5 colors

Item No.	designation	vol.	l x w x h	colour	Price
10.0390.01	Rack for 12 micro slides (without jars)	- ml	60 x 64 x 97 mm	dark gray	CHF 195.55
10.0390.02	Slide staining jar for micro slides	80 ml	76 x 64 x 92 mm	blue	CHF 100.32

10.0390.03	Slide staining jar for micro slides	80 ml	76 x 64 x 92 mm	green	CHF 100.32
10.0390.04	Slide staining jar for micro slides	80 ml	76 x 64 x 92 mm	pink	CHF 82.45
10.0390.05	Slide staining jar for micro slides	80 ml	76 x 64 x 92 mm	white	CHF 100.32
10.0390.06	Slide staining jar for micro slides	80 ml	76 x 64 x 92 mm	yellow	CHF 100.32
10.0390.07	Kit includes 5 jars (one of each color) and 1 rac	80 ml	76 x 64 x 92 mm	assorted (blue, green, pink, white, yellow)	CHF 132.00

Accessories

Item No.	Description	Price
10.0390.99	Stainless steel holder for EasyDip™	from CHF 110.62

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 07:31 oclock.

www.huberlab.ch