



## Base moulds



Product in eShop

[Home](#) / [Assortment](#) / [Research & Development](#) / [Microscopie](#) / Base moulds

### Product Data

Made of stainless steel, sizes fit most cassettes. Superior thermal exchange due to rounded corners. These moulds are for all applications in specimen embedding with Simport embedding rings and cassettes. Manufactured from high quality stainless steel for optimal thermal conductivity.

- The moulds have a highly polished surface for easy paraffin block removal
- Well corners are rounded for optimal paraffin ribboning

Item No.	dimensions (l x w x h)	Price
15.9474.01	7 x 7 x 5 mm	CHF 107.80
15.9474.02	15 x 15 x 5 mm	CHF 103.66
15.9474.03	24 x 24 x 5 mm	CHF 103.66
15.9474.04	30 x 24 x 5 mm	CHF 103.66
15.9474.05	37 x 24 x 5 mm	CHF 103.66

### More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

[www.huberlab.ch](http://www.huberlab.ch)