

Biopsy foam pads

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Microscopie](#) / [Biopsy foam pads](#)


Product Data

Made of polyester urethane foam. Biopsy foam pads are used to hold biopsies in place and prevent them from being lost during processing. They are made of a specially formulated foam which is always verified for consistency throughout in order to achieve optimum solvent flow. Biopsy samples are sandwiched between two foam pads and are placed either in tissue capsules or cassettes with metal or plastic lids.

- Can resist temperatures from 40 to +121 °C
- Blue

| Item No. | for | dimensions | Price |
|------------|----------------|--------------------|------------|
| 15.9476.01 | cassettes | 30,2 x 25,4 x 2 mm | CHF 340.11 |
| 15.9476.02 | small capsules | 25,4 x 2,7 mm | CHF 471.14 |
| 15.9476.03 | large capsules | 34,55 x 2 mm | CHF 546.11 |
| 15.9476.04 | Slimsette | 27,4 x 25,4 x 2 mm | CHF 358.96 |
| 15.9476.05 | Histosette II | 30,2 x 25,4 x 2 mm | CHF 320.57 |

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 01.05.2026 at 18:40 oclock.

www.huberlab.ch