

Stackable square base jar, with sealing cap, 5 liter, HDPE, 195x195 mm, height 215 mm,

[Product in eShop](#)[Home](#) / [Deals and news](#) / [Sale](#) /

Product Data

Weiss, stapelbar, für Crèmes, halbfeste Substanzen, Kristalle, Puder und Granulate
- Für Lebensmittel geeignet

Item No.	volume	dimensions	Price
15.1055.43	5000 ml	195x195 mm	CHF 9.00 CHF 25.30
15.1058.02	10000 ml	195x195 mm	CHF 31.50

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 02.05.2026 at 09:56 o'clock.