

Flatpack Freezer Boxen, 133x133x51mm, white

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Refrigeration and freezing](#) / [Tube boxes & freezer racks](#) /

Product Data

Wiederverwendbarer, kostengünstiger und feuchtigkeitsbeständiger Kunststoff. Diese einfach zusammensteckbaren Boxen werden in Einzelteilen geliefert um sicherzustellen, dass wertvoller Platz im Labor nicht vorzeitig belegt wird. Hergestellt aus gewelltem PP, welches langlebiger ist als Boxen aus Karton, da keine Anfälligkeit für Feuchtigkeit besteht. Wiederholtes Einfrieren und Auftauen der Boxen stellt kein Problem dar.

- Lagerung bis zu einer Temperatur von -80°C möglich
- Raster für 81 x 12 mm (max. Durchmesser) Röhrchen in einem 9 x 9 Format
- Komplettes Set mit Rastern

Item No.	dimensions (l x w x h)	Price
15.1202.62	133 x 133 x 51 mm	CHF 92.75

Accessories

Item No.	Description	Price
15.2860.13		from CHF 53.10

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 25.05.2026 at 22:31 oclock.

www.huberlab.ch