

neodisher® MediClean



Product in eShop

[Home](#) / [Assortment](#) / [Chemicals](#) / [Cleaning Chemicals](#) / [Manual cleaning](#) / neodisher® MediClean

Product Data

Mechanical and manual cleaning of medical equipment such as surgical instruments, endoscopes, anaesthetic equipment and containers.

For washer-disinfectors as well as in immersion and ultrasonic baths.

- Cleaning agent with weakly basic, cleaning supporting components in low-foaming setting.
- Residues of dried, denatured blood and blood hardened by disinfectant agents are removed by neodisher MediClean with a high degree of material protection.
- Stainless steel, instrument steel, common plastics as well as materials of optics and anaesthetic utensils are suitable for reprocessing with neodisher MediClean.
- Anodised aluminium should be checked for suitability.
- Very good cleaning effect, with perfect protection of the material, which is usually only achieved with pH-neutral products.

Item No.	execution	capacity	Price
18.9411.64	canister	10 liter	CHF 156.55
18.9411.65	canister	20 liter	CHF 279.86

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 15.05.2026 at 08:32 oclock.

www.huberlab.ch