

Mechanical multi-dispenser pipette Multipette® M4 Eppendorf

[Product in eShop](#)

Home / [Labortops_Q1_2026](#) / Mechanical multi-dispenser pipette Multipette® M4 Eppendorf

Product Data

Say goodbye to tedious, time-consuming serial pipetting! The Multipette M4 pipette from Eppendorf is ideal for repetitive dispensing tasks in long series, eliminating the need for repeated liquid aspiration. Paired with Combitips® advanced dispenser tips, the Multipette effortlessly handles viscous, volatile or foaming liquids that pose challenges to traditional pipettes. The hermetically sealed tips in the Multipette/Combitips positive displacement pipetting system ensure a high level of safety when working with infectious, toxic, or radioactive liquids.

- Mechanical positive displacement multi-dispenser: fill the tip once and dispense up to 100 times without a refill
- Automatic Combitips® advanced dispenser tip recognition eliminates time-consuming volume calculations
- Easy-to-read display guarantees stress-free and intuitive operation with a digital display of the dispensing volume
- Wide dispensing range: 1 µL to 10 mL for flexible operation
- Stress-free work via integrated step counter: dispensing procedures can be continued error free after an interruption or distraction
- Integrated »Sleep« function switches the instrument off, when it is not being used, thus reducing energy consumption
- Fully emptied Combitips® advanced tip can be easily ejected with one hand using the operating lever
- Positive displacement principle: Ideal for viscous or foaming solutions and liquids with high vapor pressure
- Safe handling of toxic, radioactive or infectious liquids
- Eco-conscious production: Multipette M4 and Combitips® advanced tips are ACT®-labeled to help you make an informed choice when it comes to sustainability

Item No.	designation	Price
11.4988.09	Multipette® M4, 1-channel, 1 µL – 10 mL, incl. 1x Combitips® advanced 2.5 mL (green) tip and Pipette Holder 2	CHF 506.00

Accessories

Item No.	Description	Price
11.4807.06	Accessories for pipettes Eppendorf	from CHF 28.40
11.4990.20	Combitips® advanced Eppendorf	from CHF 18.85
11.3116.31	Pipette Holder System Eppendorf	from CHF 40.70

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 25.04.2026 at 21:19 o'clock.

www.huberlab.ch