

EDISONITE® CLASSIC Universal detergents

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Safety & Hygiene](#) / [Cleaning](#) / [Cleansers & soaps](#) / EDISONITE® CLASSIC Universal detergents

Product Data

Alkaline detergent. Powder compound for immersion and ultrasonic bath. Ideal for removal of natural fats and oils, proteins and blood. Durable protection against corrosion damage.

- **Ingredients:** > 30% phosphates, < 5% anionic surfactants, non-ionic surfactants, corrosion inhibitors, auxiliary agents
- **Concentration for use:** 0.5 - 1% (5 - 10 g/l).
- **Chemical and physical specifications:** Lightgreen powder
- **Bulk density:** approx. 850 kg/m³.
- **pH value** (0.5%): approx. 11.5. pH value (1%): approx. 11.5
- Use on zinc and non-anodized aluminium is restricted

Item No.	capacity	Price
18.1101.60	5 kg bucket	CHF 234.35
18.1101.61	1 kg jar	CHF 58.10

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 08.05.2026 at 14:12 o'clock.