

# Pursept® A Xpress Surface Disinfecting Detergents

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Safety & Hygiene](#) / [Disinfection](#) / Pursept® A Xpress Surface Disinfecting Detergents

## Product Data

Refillable pump-atomizer without gas. Pleasant, fresh fragrance. Effective in only 15 seconds! DGHM/VAH tested and RKL compliant, with aerosol-free application (disinfectant wipes).

- **Ingredients:** 100 g Pursept® A Xpress contain: 55 g ethanol, 0.03 g n-alkyl-aminopropyl-glycine, auxiliary agents
- **Chemical and physical specifications:**  
Clear, colorless liquid  
Density: approx. 0.9 g/ml  
pH value: approx. 7,5
- **Mikrobiologische Wirksamkeit/Anwendungskonzentration:**  
Flächendesinfektion:  
Bakterien und Pilze bei hoher organischer Belastung gem. DGHM/VAH: 0,25% (2,5 g/l)/4 Std., 0,5% (5 g/l)/1 Std., 1% (10 g/l)/30 Min., 2% (20 g/l)/15 Min.  
Begrenzt viruzid\* (inkl. HIV, HBV, HCV): 1% (10 g/l)/15 Min.  
Inaktivierung humane Rotaviren: 0,25% (2,5 g/l)/2 Min.  
Wirksamkeit gegenüber TbB: 2% (20 g/l)/1 Std.

Item No.	capacity	Price
18.4487.03	atomizer head	CHF 35.40

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 12.05.2026 at 13:57 o'clock.