

## Pursept® Wipes XL Fleece wipes



Product in eShop

[Home](#) / [Assortment](#) / [Safety & Hygiene](#) / [Disinfection](#) / Pursept® Wipes XL Fleece wipes

### Product Data

Fleece wipes for impregnation with universal or disinfecting detergents. The lint-free and tear-proof wipes are designed for both optimal disinfectant delivery and excellent contaminant uptake. Due to their size they are particularly suitable to disinfect your workplace, lab table or clean bench quickly and reliably. The wipes can be stored in the convenient dispenser box. For high-volume use, wipes can be soaked with liquid disinfectant directly in the handy dispenser box (lasts up to 28 days). We recommend the usage of the surface disinfecting detergent Pursept® A Xpress.

#### - Technical Information:

- Fleece wipes made of 100% polyester, dry
- 70 wipes (230 x 280 mm) per roll
- Dispenser Box, refillable and firmly closable

Item No.	capacity	Price
18.4487.01	Dispenser box	CHF 35.90

### More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 03.05.2026 at 00:12 o'clock.

[www.huberlab.ch](http://www.huberlab.ch)