


Contact dishes

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Life Science](#) / [Microbiology](#) / [Petri dishes](#) / [Contact dishes](#)

Product Data

Made of PS, transparent. Sterile (gamma-sterilised). For detecting microorganisms and testing the effectiveness of cleaning and disinfection on flat surfaces using contact cultures.

- Dimensions: ? 65 x H 15 mm
- Optimised dish design produces a meniscus on the surface of the agar when the contact dish is filled
- Optionally available with vents for improved gas exchange, or without vents
- Graduated and divided into 10 mm squares to allow quick and easy analysis

Item No.	vol.	vents	packaging	Price
7.9210.17	20 ml	-	25 packs x 20 pcs	CHF 251.08
7.9210.18	15 ml	+	25 packs x 20 pcs	CHF 412.50

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 10.05.2026 at 10:33 o'clock.