

Manual colony counter Scan® 100 Interscience

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Life Science](#) / [Cell Biology](#) / [Colony counters](#) / [Manual colony counter Scan® 100 Interscience](#)

Product Data

Scan® 100 is a high-tech manual colony counter remarkable for its user-friendly functions and wonderful ergonomics. A LED lighting system with Dark-Field technology gives you an accurate and contrasted view of the colonies. The integrated USB port allows export of the results to guarantee traceability.

- For all Petri dishes: diameter 55, 65, 90, 150 mm - Multi-media
- Dimensions (L x W x H): 31 x 25 x 18 cm
- Weight: 4.5 kg
- Stainless steel
- Adjustable sensitivity of touch surface
- Multifunction digital screen
- Counting: from 0 to 1999 CFU
- USB port: easy data transfer via any software (PC/Mac)
- Dark Field technology for a contrasted view of the colonies
- Adjustable counting beeps: 4 sound levels
- Adjustable LED intensity
- Adapted to any pen
- Built-in storage rack

Item No.	description	Price
17.6072.10	Scan® 100	CHF 1'987.00
17.6072.15	Magnifying lens	CHF 407.70

Accessories

Item No.	Description	Price
17.4370.10	Printer dataLink® pro interscience	from CHF 11'882.20

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 22:00 o'clock.

www.huberlab.ch